

Product Name: ZXDC Mouse Monoclonal Antibody
Catalog #: AMM83346



Summary

Production Name	ZXDC Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ZXDC
Alternative Names	ZXDC
Gene ID	
SwissProt ID	Q2QGD7. Purified recombinant fragment of human ZXDC

Application

Dilution Ratio	IHC:1:500
Molecular Weight	/


Background

Product Name: ZXDC Mouse Monoclonal Antibody
Catalog #: AMM83346



Research Area

Image Data



IHC-P analysis of cervical cancer tissue by ZXDC antibody (AMM83346)

Note

For research use only.