
Product Name: PFKM Mouse Monoclonal Antibody**Catalog #: AMM83349**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	PFKM
Alternative Names	PFKM
Gene ID	
SwissProt ID	P08237
Immunogen	Purified recombinant fragment of human PFKM

Background**Research Area**

Image Data



Lane 1: SH-SY5Y lysate ; Lane 2: PFKM immunoprecipitated from SH-SY5Y lysate by PFKM Mouse Monoclonal Antibody ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: PFKM Mouse Monoclonal Antibody can immunoprecipitate PFKM;



IHC-P analysis of skeletal muscle tissue by PFKM antibody (AMM83349)