

Product Name: KIFBP Mouse Monoclonal Antibody
Catalog #: AMM83360



Summary

Production Name	KIFBP Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IP,IHC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	KIFBP
Alternative Names	KIFBP
Gene ID	
SwissProt ID	Q96EK5.Purified recombinant fragment of human KIFBP

Application

Dilution Ratio	IP:1:100-1:200,IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: KIFBP Mouse Monoclonal Antibody
Catalog #: AMM83360




Research Area

Image Data



Lane 1: Mouse kidney tissue lysate ; Lane 2: KIFBP immunoprecipitated from mouse kidney tissue lysate by AMM83360 ;
Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83360 can immunoprecipitate KIFBP;



IHC-P analysis of kidney tissue by KIFBP antibody (AMM83360)

Note

For research use only.