

Summary

Production Name	BPGM Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	BPGM
Alternative Names	BPGM
Gene ID	
SwissProt ID	P07738. Purified recombinant fragment of human BPGM

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/

Background

Product Name: BPGM Mouse Monoclonal Antibody
Catalog #: AMM83366



Research Area

Image Data



IHC-P analysis of placenta tissue by BPGM antibody (AMM83366)

Note

For research use only.