
Product Name: FGF12 Mouse Monoclonal Antibody**Catalog #: AMM83372**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	27kDa

Antigen Information


Gene Name	FGF12
Alternative Names	FGF12
Gene ID	
SwissProt ID	P61328
Immunogen	Purified recombinant fragment of human FGF12

Background**Research Area**

Image Data



Lane 1: 15 µg of rat brain tissue lysate; Lane 2: 15 µg of mouse brain tissue lysate; Lane 3: 15 µg of brain tissue lysate; Result: FGF12 Mouse Monoclonal Antibody can detect FGF12 by Western blotting;



Lane 1: Rat brain lysate ; Lane 2: FGF12 immunoprecipitated from rat brain lysate by AMM83372 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ; Result: AMM83372 can immunoprecipitate FGF12;