
Product Name: TP53I3 Mouse Monoclonal Antibody**Catalog #: AMM83376**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	TP53I3
Alternative Names	TP53I3
Gene ID	
SwissProt ID	Q53FA7
Immunogen	Purified recombinant fragment of human TP53I3

Background**Research Area**

Image Data



Lane 1: Hep G2 lysate ; Lane 2: TP53I3 immunoprecipitated from Hep G2 lysate by TP53I3 Mouse Monoclonal Antibody ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: TP53I3 Mouse Monoclonal Antibody can immunoprecipitate TP53I3;