
Product Name: TSC1 Mouse Monoclonal Antibody**Catalog #: AMM83377**

For research use only.

Summary

| | |
|----------------------|---|
| Description | Mouse monoclonal Antibody |
| Host | Mouse |
| Application | IP |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Form | Liquid |
| Concentration | 1mg/ml |
| Storage | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| Shipping | Ice bags |
| Buffer | Purified antibody in PBS with 0.05% sodium azide |
| Purification | Affinity Purification |

Application


| | |
|-------------------------|----------------|
| Dilution Ratio | IP 1:100-1:200 |
| Molecular Weight | / |

Antigen Information

| | |
|--------------------------|---|
| Gene Name | TSC1 |
| Alternative Names | TSC1 |
| Gene ID | |
| SwissProt ID | Q92574 |
| Immunogen | Purified recombinant fragment of human TSC1 |

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate ; Lane 2: TSC1 immunoprecipitated from HEK-293 lysate by TSC1 Mouse Monoclonal Antibody ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: TSC1 Mouse Monoclonal Antibody can immunoprecipitate TSC1;