

Product Name: HERC5 Mouse Monoclonal Antibody**Catalog #: AMM83385**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	116.9kDa

Antigen Information


Gene Name	HERC5
Alternative Names	HERC5
Gene ID	
SwissProt ID	Q9UII4
Immunogen	Purified recombinant fragment of human HERC5

Background**Research Area**


Image Data



Western blot analysis using HERC5 mouse mAb against HEK-293 lysate



Lane 1: HEK-293 lysate ; Lane 2: HERC5 immunoprecipitated from HEK-293 lysate by AMM83385 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83385 can immunoprecipitate HERC5;



IHC-P analysis of testis tissue by HERC5 antibody (AMM83385)