
Product Name: ARID1A Mouse Monoclonal Antibody**Catalog #: AMM83386**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	242kDa

Antigen Information


Gene Name	ARID1A
Alternative Names	ARID1A
Gene ID	
SwissProt ID	O14497
Immunogen	Purified recombinant fragment of human ARID1A

Background**Research Area**

Image Data



Western blot analysis using ARID1A mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: ARID1A immunoprecipitated from Jurkat lysate by AMM83386 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83386 can immunoprecipitate ARID1A;