
Product Name: HERC4 Mouse Monoclonal Antibody**Catalog #: AMM83389**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	118.6kDa

Antigen Information

Gene Name	HERC4
Alternative Names	HERC4
Gene ID	
SwissProt ID	Q5GLZ8
Immunogen	Purified recombinant fragment of human HERC4

Background**Research Area**

Image Data



Western blot analysis using HERC4 mouse mAb against rat liver tissue lysate