
Product Name: IL25 Mouse Monoclonal Antibody**Catalog #: AMM83398**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	20.3kDa

Antigen Information

Gene Name	IL25
Alternative Names	IL25
Gene ID	
SwissProt ID	Q9H293
Immunogen	Purified recombinant fragment of human IL25

Background**Research Area**

Image Data

