
Product Name: CA3 Mouse Monoclonal Antibody**Catalog #: AMM83402**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	38kDa

Antigen Information

Gene Name	CA3
Alternative Names	CA3
Gene ID	
SwissProt ID	P07451
Immunogen	Purified recombinant fragment of human CA3

Background**Research Area**

Image Data



Western blot analysis using CA3 mouse mAb against mouse skeletal muscle tissue lysate