

Product Name: TRIM72 Mouse Monoclonal Antibody
Catalog #: AMM83404



Summary

Production Name	TRIM72 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	TRIM72
Alternative Names	TRIM72
Gene ID	
SwissProt ID	Q6ZMU5. Purified recombinant fragment of human TRIM72

Application

Dilution Ratio	IP:1:100-1:200
Molecular Weight	/


Background

Product Name: TRIM72 Mouse Monoclonal Antibody
Catalog #: AMM83404



Research Area

Image Data



Lane 1: Rat skeletal muscle tissue lysate; Lane 2: TRIM72 immunoprecipitated from rat skeletal muscle tissue lysate by AMM83404; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83404 can immunoprecipitate TRIM72;

Note

For research use only.