
Product Name: LONRF2 Mouse Monoclonal Antibody**Catalog #: AMM83410**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	LONRF2
Alternative Names	LONRF2
Gene ID	
SwissProt ID	Q1L5Z9
Immunogen	Purified recombinant fragment of human LONRF2

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate; Lane 2: LONRF2 immunoprecipitated from HEK-293 lysate by LONRF2 Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: LONRF2 Mouse Monoclonal Antibody can immunoprecipitate LONRF2;