
Product Name: TRAIP Mouse Monoclonal Antibody**Catalog #: AMM83411**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	TRAIP
Alternative Names	TRAIP
Gene ID	
SwissProt ID	Q9BWF2
Immunogen	Purified recombinant fragment of human TRAIP

Background**Research Area**

Image Data



Lane 1: Mouse brain tissue lysate; Lane 2: TRAIP immunoprecipitated from mouse brain tissue lysate by TRAIP Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: TRAIP Mouse Monoclonal Antibody can immunoprecipitate TRAIP;