

Product Name: AREL1 Mouse Monoclonal Antibody**Catalog #: AMM83417**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:25-1:200
Molecular Weight	/

Antigen Information

Gene Name	AREL1
Alternative Names	AREL1
Gene ID	
SwissProt ID	O15033
Immunogen	Purified recombinant fragment of human AREL1

Background**Research Area**

Image Data



IHC-P analysis of tonsil tissue by AREL1 antibody (AREL1 Mouse Monoclonal Antibody)