

Product Name: CXCL13 Mouse Monoclonal Antibody
Catalog #: AMM83418



Summary

Production Name	CXCL13 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CXCL13
Alternative Names	CXCL13
Gene ID	
SwissProt ID	O43927. Purified recombinant fragment of human CXCL13

Application

Dilution Ratio	IHC:1:20-1:100
Molecular Weight	/


Background

Product Name: CXCL13 Mouse Monoclonal Antibody
Catalog #: AMM83418



Research Area

Image Data



IHC-P analysis of tonsil tissue by CXCL13 antibody (AMM83418)

Note

For research use only.