
Product Name: Rapgef4 Mouse Monoclonal Antibody**Catalog #: AMM83430**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:25-1:200
Molecular Weight	/

Antigen Information

Gene Name	Rapgef4
Alternative Names	Rapgef4
Gene ID	
SwissProt ID	Q9EQZ6
Immunogen	Purified recombinant fragment of mouse Rapgef4

Background**Research Area**

Image Data



IHC-P analysis of mouse liver tissue by Rapgef4 antibody (Rapgef4 Mouse Monoclonal Antibody)