

Product Name: VASH1 Mouse Monoclonal Antibody
Catalog #: AMM83443



Summary

Production Name	VASH1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	VASH1
Alternative Names	VASH1
Gene ID	
SwissProt ID	Q7L8A9. Purified recombinant fragment of human VASH1

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	40kDa


Background

Product Name: VASH1 Mouse Monoclonal Antibody
Catalog #: AMM83443



Research Area

Image Data



Western blot analysis using VASH1 mouse mAb against Jurkat lysate and 15 μ g

Note

For research use only.