

Product Name: NDUF8 Mouse Monoclonal Antibody**Catalog #: AMM83453**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:25-1:200,IP 1:100-1:200
Molecular Weight	21.8kDa

Antigen Information

Gene Name	NDUF8
Alternative Names	NDUF8
Gene ID	
SwissProt ID	O95169
Immunogen	Purified recombinant fragment of human NDUF8

Background**Research Area**

Image Data



Lane 1: 15 µg of HeLa lysate