
Product Name: RECQL Mouse Monoclonal Antibody**Catalog #: AMM83466**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:20-1:100,IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	RECQL
Alternative Names	RECQL
Gene ID	
SwissProt ID	P46063
Immunogen	Purified recombinant fragment of human RECQL

Background**Research Area**

Image Data



IHC-P analysis of heart muscle tissue by RECQL antibody (RECQL Mouse Monoclonal Antibody)