

Product Name: CDH15 Mouse Monoclonal Antibody**Catalog #: AMM83473**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	89kDa

Antigen Information

Gene Name	CDH15
Alternative Names	CDH15
Gene ID	
SwissProt ID	P55291
Immunogen	Purified recombinant fragment of human CDH15

Background**Research Area**

Image Data



Western blot analysis using CDH15 mouse mAb against HeLa lysate and mouse brain tissue lysate