

Product Name: HPCAL1 Mouse Monoclonal Antibody**Catalog #: AMM83475**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	22.3kDa

Antigen Information

Gene Name	HPCAL1
Alternative Names	HPCAL1
Gene ID	
SwissProt ID	P37235
Immunogen	Purified recombinant fragment of human HPCAL1

Background**Research Area**

Image Data



Western blot analysis using HPCAL1 mouse mAb against rat brain tissue lysate and mouse brain tissue lysate