

Product Name: ZXDC Mouse Monoclonal Antibody
Catalog #: AMM83477



Summary

Production Name	ZXDC Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ZXDC
Alternative Names	ZXDC
Gene ID	
SwissProt ID	Q2QGD7.Purified recombinant fragment of human ZXDC

Application

Dilution Ratio	WB:1:2000,IP:1:100-1:200
Molecular Weight	90kDa


Background

Product Name: ZXDC Mouse Monoclonal Antibody
Catalog #: AMM83477




Research Area

Image Data



Western blot analysis using ZXDC mouse mAb against mouse kidney tissue lysate



Lane 1: Mouse kidney tissue lysate; Lane 2: ZXDC immunoprecipitated from mouse kidney tissue by AMM83477; Lane 3:
The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83477 can immunoprecipitate ZXDC;

Note

For research use only.