

Product Name: HSPA2 Mouse Monoclonal Antibody**Catalog #: AMM83480**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	HSPA2
Alternative Names	HSPA2
Gene ID	
SwissProt ID	P54652
Immunogen	Purified recombinant fragment of human HSPA2

Background**Research Area**

Image Data



Lane 1: Mouse testis tissue lysate; Lane 2: HSPA2 immunoprecipitated from mouse testis tissue lysate by HSPA2 Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: HSPA2 Mouse Monoclonal Antibody can immunoprecipitate HSPA2;