
Product Name: FBXO34 Mouse Monoclonal Antibody**Catalog #: AMM83495**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:20-1:100
Molecular Weight	/

Antigen Information

Gene Name	FBXO34
Alternative Names	FBXO34
Gene ID	
SwissProt ID	Q9NWN3
Immunogen	Purified recombinant fragment of human FBXO34

Background**Research Area**

Image Data



IHC-P analysis of thyroid gland cancer tissue by FBXO34 antibody (FBXO34 Mouse Monoclonal Antibody)