
Product Name: CD2AP Mouse Monoclonal Antibody**Catalog #: AMM83502**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	CD2AP
Alternative Names	CD2AP
Gene ID	
SwissProt ID	Q9Y5K6
Immunogen	Purified recombinant fragment of human CD2AP

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate; Lane 2: CD2AP immunoprecipitated from HEK-293 lysate by CD2AP Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: CD2AP Mouse Monoclonal Antibody can immunoprecipitate CD2AP;