
Product Name: PSMA1 Mouse Monoclonal Antibody**Catalog #: AMM83503**

For research use only.

Summary

| | |
|----------------------|---|
| Description | Mouse monoclonal Antibody |
| Host | Mouse |
| Application | WB |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | Mouse IgG2b |
| Clonality | Monoclonal |
| Form | Liquid |
| Concentration | 1mg/ml |
| Storage | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| Shipping | Ice bags |
| Buffer | Purified antibody in PBS with 0.05% sodium azide |
| Purification | Affinity Purification |

Application


| | |
|-------------------------|------------------|
| Dilution Ratio | WB 1:1000-1:2000 |
| Molecular Weight | 29.6kDa |

Antigen Information

| | |
|--------------------------|--|
| Gene Name | PSMA1 |
| Alternative Names | PSMA1 |
| Gene ID | |
| SwissProt ID | P25786 |
| Immunogen | Purified recombinant fragment of human PSMA1 |

Background**Research Area**

Image Data



Western blot analysis using PSMA1 mouse mAb against HeLa lysate and 15 μ g