

**Product Name: PPP1R27 Mouse Monoclonal Antibody**  
**Catalog #: AMM83504**



## Summary

<b>Production Name</b>	PPP1R27 Mouse Monoclonal Antibody
<b>Description</b>	Mouse Monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	IHC
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	PPP1R27
<b>Alternative Names</b>	PPP1R27
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q86WC6. Purified recombinant fragment of human PPP1R27

## Application

<b>Dilution Ratio</b>	IHC:1:25-1:200
<b>Molecular Weight</b>	/

## Background


**Product Name: PPP1R27 Mouse Monoclonal Antibody**  
**Catalog #: AMM83504**



---

## Research Area

## Image Data



IHC-P analysis of skeletal muscle tissue by PPP1R27 antibody (AMM83504)

## Note

For research use only.