

Product Name: MCF2L2 Mouse Monoclonal Antibody
Catalog #: AMM83506



Summary

Production Name	MCF2L2 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	MCF2L2
Alternative Names	MCF2L2
Gene ID	
SwissProt ID	Q86YR7. Purified recombinant fragment of human MCF2L2

Application

Dilution Ratio	IHC:1:20-1:100
Molecular Weight	/


Background

Product Name: MCF2L2 Mouse Monoclonal Antibody
Catalog #: AMM83506



Research Area

Image Data



IHC-P analysis of brain tissue by MCF2L2 antibody (AMM83506)

Note

For research use only.