

Product Name: PDP1 Mouse Monoclonal Antibody
Catalog #: AMM83507



Summary

Production Name	PDP1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	PDP1
Alternative Names	PDP1
Gene ID	
SwissProt ID	Q9P0J1. Purified recombinant fragment of human PDP1

Application

Dilution Ratio	IP:1:100-1:200
Molecular Weight	/


Background

Product Name: PDP1 Mouse Monoclonal Antibody
Catalog #: AMM83507



Research Area

Image Data



Lane 1: Jurkat lysate; Lane 2: PDP1 immunoprecipitated from Jurkat lysate by AMM83507; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83507 can immunoprecipitate PDP1;

Note

For research use only.