

Product Name: ELAVL4 Mouse Monoclonal Antibody
Catalog #: AMM83508



Summary

Production Name	ELAVL4 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ELAVL4
Alternative Names	ELAVL4
Gene ID	
SwissProt ID	P26378-6. Purified recombinant fragment of human ELAVL4

Application

Dilution Ratio	IHC:1:20-1:100
Molecular Weight	/


Background

Product Name: ELAVL4 Mouse Monoclonal Antibody
Catalog #: AMM83508



Research Area

Image Data



IHC-P analysis of cerebral cortex tissue by ELAVL4 antibody (AMM83508)

Note

For research use only.