

Product Name: VWA5A Mouse Monoclonal Antibody
Catalog #: AMM83509



Summary

Production Name	VWA5A Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	VWA5A
Alternative Names	VWA5A
Gene ID	
SwissProt ID	O00534. Purified recombinant fragment of human VWA5A

Application

Dilution Ratio	IHC:1:20-1:100
Molecular Weight	/


Background

Product Name: VWA5A Mouse Monoclonal Antibody
Catalog #: AMM83509



Research Area

Image Data



IHC-P analysis of heart muscle tissue by VWA5A antibody (AMM83509)

Note

For research use only.