

Product Name: CLEC1A Mouse Monoclonal Antibody
Catalog #: AMM83511



Summary

Production Name	CLEC1A Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CLEC1A
Alternative Names	CLEC1A
Gene ID	
SwissProt ID	Q8NC01. Purified recombinant fragment of human CLEC1A

Application

Dilution Ratio	IHC:1:25-1:200
Molecular Weight	/


Background

Product Name: CLEC1A Mouse Monoclonal Antibody
Catalog #: AMM83511



Research Area

Image Data



IHC-P analysis of lung tissue by CLEC1A antibody (AMM83511)

Note

For research use only.