

Product Name: IFIT3 Mouse Monoclonal Antibody
Catalog #: AMM83512



Summary

Production Name	IFIT3 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	IFIT3
Alternative Names	IFIT3
Gene ID	
SwissProt ID	O14879. Purified recombinant fragment of human IFIT3

Application

Dilution Ratio	WB:1:2000
Molecular Weight	56kDa


Background

Product Name: IFIT3 Mouse Monoclonal Antibody
Catalog #: AMM83512



Research Area

Image Data



Western blot analysis using IFIT3 mouse mAb against HT-29 lysate

Note

For research use only.