
Product Name: PIK3AP1 Mouse Monoclonal Antibody**Catalog #: AMM83518**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	PIK3AP1
Alternative Names	PIK3AP1
Gene ID	
SwissProt ID	Q6ZUJ8
Immunogen	Purified recombinant fragment of human PIK3AP1

Background**Research Area**

Image Data



IHC-P analysis of tonsil tissue by PIK3AP1 antibody (PIK3AP1 Mouse Monoclonal Antibody)