
Product Name: NPTXR Mouse Monoclonal Antibody**Catalog #: AMM83527**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	NPTXR
Alternative Names	NPTXR
Gene ID	
SwissProt ID	O95502
Immunogen	Purified recombinant fragment of human NPTXR

Background**Research Area**

Image Data



Lane 1: Mouse brain tissue lysate; Lane 2: NPTXR immunoprecipitated from mouse brain tissue lysate by NPTXR Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: NPTXR Mouse Monoclonal Antibody can immunoprecipitate NPTXR;