
Product Name: PTN Mouse Monoclonal Antibody**Catalog #: AMM83530**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	PTN
Alternative Names	PTN
Gene ID	
SwissProt ID	P21246
Immunogen	Purified recombinant fragment of human PTN

Background**Research Area**

Image Data



Lane 1: NIH/3T3 lysate; Lane 2: PTN immunoprecipitated from NIH/3T3 lysate by PTN Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: PTN Mouse Monoclonal Antibody can immunoprecipitate PTN;