
Product Name: TRIM43B Mouse Monoclonal Antibody**Catalog #: AMM83532**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	TRIM43B
Alternative Names	TRIM43B
Gene ID	
SwissProt ID	A6NCK2
Immunogen	Purified recombinant fragment of human TRIM43B

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate (transfected with DYKDDDDK tagged TRIM43B (Cat. P69024QATN) at N-terminus); Lane 2: TRIM43B immunoprecipitated from the same lysate as Lane 1; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: TRIM43B Mouse Monoclonal Antibody can immunoprecipitate TRIM43B;