
Product Name: RNF113B Mouse Monoclonal Antibody**Catalog #: AMM83533**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	RNF113B
Alternative Names	RNF113B
Gene ID	
SwissProt ID	Q8IZP6
Immunogen	Purified recombinant fragment of human RNF113B

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate (transfected with DYKDDDDK tagged RNF113B (Cat. P65017QATN) at N-terminus); Lane 2: TRIM49 immunoprecipitated from the same lysate as Lane 1; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: RNF113B Mouse Monoclonal Antibody can immunoprecipitate RNF113B;



Lane 1: HEK-293 lysate (transfected with EGFP tagged RNF113B (Cat. P65017QBTN) at C-terminus); Lane 2: TRIM49 immunoprecipitated from the same lysate Lane 1; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83533 can immunoprecipitate RNF113B;