

Product Name: PSMC3IP Mouse Monoclonal Antibody**Catalog #: AMM83538**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:20-1:100
Molecular Weight	/

Antigen Information

Gene Name	PSMC3IP
Alternative Names	PSMC3IP
Gene ID	
SwissProt ID	Q9P2W1
Immunogen	Purified recombinant fragment of human PSMC3IP

Background**Research Area**

Image Data



IHC-P analysis of testis tissue by PSMC3IP antibody (PSMC3IP Mouse Monoclonal Antibody)