
Product Name: PRKAR1B Mouse Monoclonal Antibody**Catalog #: AMM83540**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:50-1:200
Molecular Weight	/

Antigen Information

Gene Name	PRKAR1B
Alternative Names	PRKAR1B
Gene ID	
SwissProt ID	P31321
Immunogen	Purified recombinant fragment of human PRKAR1B

Background**Research Area**

Image Data



IHC-P analysis of brain tissue by PRKAR1B antibody (PRKAR1B Mouse Monoclonal Antibody)