

Product Name: CXCL6 Mouse Monoclonal Antibody**Catalog #: AMM83552**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	11.9kDa

Antigen Information

Gene Name	CXCL6
Alternative Names	CXCL6
Gene ID	
SwissProt ID	P80162
Immunogen	Purified recombinant fragment of human CXCL6

Background**Research Area**

Image Data



Result: CXCL6 Mouse Monoclonal Antibody can detect CXCL6 by Western blotting.