

Product Name: SPINT1 Mouse Monoclonal Antibody
Catalog #: AMM83553



Summary

Production Name	SPINT1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	SPINT1
Alternative Names	SPINT1
Gene ID	
SwissProt ID	O43278. Purified recombinant fragment of human SPINT1

Application

Dilution Ratio	IHC:1:25-1:200
Molecular Weight	/


Background

Product Name: SPINT1 Mouse Monoclonal Antibody
Catalog #: AMM83553



Research Area

Image Data



IHC-P analysis of esophagus tissue by SPINT1 antibody (AMM83553)

Note

For research use only.