

---

**Product Name: TNRC18 Mouse Monoclonal Antibody****Catalog #: AMM83558**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	IHC
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**


<b>Dilution Ratio</b>	IHC 1:50-1:100
<b>Molecular Weight</b>	/

**Antigen Information**

<b>Gene Name</b>	TNRC18
<b>Alternative Names</b>	TNRC18
<b>Gene ID</b>	
<b>SwissProt ID</b>	O15417
<b>Immunogen</b>	Purified recombinant fragment of human TNRC18

**Background****Research Area**

## Image Data



IHC-P analysis of placenta tissue by TNRC18 antibody (TNRC18 Mouse Monoclonal Antibody)