
Product Name: AARSD1 Mouse Monoclonal Antibody**Catalog #: AMM83567**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:20-1:100
Molecular Weight	/

Antigen Information

Gene Name	AARSD1
Alternative Names	AARSD1
Gene ID	
SwissProt ID	Q9BTE6
Immunogen	Purified recombinant fragment of human AARSD1

Background**Research Area**

Image Data



IHC-P analysis of testis tissue by AARSD1 antibody (AARSD1 Mouse Monoclonal Antibody)