

Product Name: RIPK3 Mouse Monoclonal Antibody
Catalog #: AMM83568



Summary

Production Name	RIPK3 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	RIPK3
Alternative Names	RIPK3
Gene ID	
SwissProt ID	Q9Y572. Purified recombinant fragment of human RIPK3

Application

Dilution Ratio	IHC:1:25-1:200
Molecular Weight	/


Background

Product Name: RIPK3 Mouse Monoclonal Antibody
Catalog #: AMM83568



Research Area

Image Data



IHC-P analysis of pancreas tissue by RIPK3 antibody (AMM83568)

Note

For research use only.