
Product Name: MYH9 Mouse Monoclonal Antibody**Catalog #: AMM83577**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:250,IP 1:100-1:200
Molecular Weight	227kDa

Antigen Information


Gene Name	MYH9
Alternative Names	MYH9
Gene ID	
SwissProt ID	P35579
Immunogen	Purified recombinant fragment of human MYH9

Background**Research Area**


Image Data



Lane 1: 15 μ g of HeLa lysate; Lane 2: 15 μ g of HT-29 lysate; Lane 3: 15 μ g of A-431 lysate; Result: MYH9 Mouse Monoclonal Antibody can detect MYH9 by Western blotting;



Lane 1: HeLa lysate; Lane 2: MYH9 immunoprecipitated from HeLa lysate by AMM83577; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83577 can immunoprecipitate MYH9;



IHC-P analysis of kidney tissue by MYH9 antibody (AMM83577)