

Product Name: Human Lambda Mouse Monoclonal Antibody
Catalog #: AMM83585

Summary

Production Name	Human Lambda Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Human Lambda
Alternative Names	Human Lambda
Gene ID	
SwissProt ID	POCG04. Human Lambda

Application

Dilution Ratio	IHC:1:500
Molecular Weight	/


Background

Product Name: Human Lambda Mouse Monoclonal Antibody
Catalog #: AMM83585



Research Area

Image Data



IHC-P analysis of tonsil tissue by Lambda antibody (AMM83585)

Note

For research use only.