
Product Name: INHBA Mouse Monoclonal Antibody**Catalog #: AMM83589**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgA
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:250
Molecular Weight	52kDa

Antigen Information


Gene Name	INHBA
Alternative Names	INHBA
Gene ID	
SwissProt ID	P08476
Immunogen	Purified recombinant fragment of human INHBA

Background**Research Area**

Image Data



Result: INHBA Mouse Monoclonal Antibody can detect INHBA by Western blotting.



IHC-P analysis of lung cancer tissue by INHBA antibody (AMM83589)